

Glucose Vision

Blood Glucose Monitoring System



NEW!

★ **Accurate**

★ **Easy to Read**

★ **Small and portable**

★ **Large memory storage**

★ **Only 1 μ l blood sample required**

**Results in just
5 seconds**

Glucose Vision

Blood Glucose Monitoring System

Kit:

1. Glucose Meter
2. Test Strips
3. Adjustable Lancet
4. Positive Control Reagent




Procedure:



1. Insert Test Strip



2. Apply Blood to Test Spot



3. Read Results in 5 seconds

Specifications:

- ☆ **Test Time:**
5 seconds
- ☆ **Blood Volume:**
Only 1 micro liter
- ☆ **Size:**
3.2" (L) x 2.3" (W) x 0.8 (H)
8.1cm (L) x 5.8cm (W) x 2.0cm (H)
- ☆ **Weight:**
2.54 ounces (72g) including 2 batteries
- ☆ **Memory Capacity:**
Up to 300 recent results
- ☆ **Calibration Curve:**
54 total code numbers
- ☆ **Blood Source:**
Fresh whole blood
- ☆ **Display:**
LCD Display
- ☆ **System Operating Range:**
Temperature: 10 - 35 C (50 - 95 F)
Lower results are displayed as "tLo"
Higher results are displayed as "tLi"
Relative Humidity: 25 - 90%
- ☆ **Measurement Range:**
20 to 600 mg/dL (1.1 to 33.3 mMoie/L)
Lower results are displayed as "LO"
Higher results are displayed as "HI"
- ☆ **Power Supply:**
2 replaceable AAA size batteries



2774 Loker Av. West
Calrsbad - CA, 92108

Phone: (760) 930-2200
Fax: (760) 930-2212

Email: marketing@syntron.net
Website: www.syntron.net